



Kirsch Sinus Lift Instrument Set SMALL

for sinus floor elevation, modif. acc. to Dr. Kai Zwanzig





24.995.00S

Kirsch Sinus Lift Instrument Set SMALL,
modif. acc. to Dr. Kai Zwanzig,
incl. 1 x washbasket with lid 85.195.00 and
arrangement for Sinus Lift Instrument Set

Complete Set consisting of:



85.251.04

Mixing Cup, stainless steel,
with plastic lid, Ø 40 mm

85.251.14

Mixing Cup without plastic lid,
stainless steel, Ø 40 mm



41.868.01S

Kirsch, Sinus 1, SMALL, Elevator 25°,
Titanium **ZEPF**-Line Handle, double-ended



41.868.02S

Kirsch, Sinus 2, SMALL, Elevator 60°/60°,
Titanium **ZEPF**-Line Handle, double-ended



41.868.03S

Kirsch, Sinus 3, SMALL, Elevator 90°,
Titanium **ZEPF**-Line Handle, double-ended



41.868.04S

Kirsch, Sinus 4, SMALL, Elevator 120°/120°,
Titanium **ZEPF**-Line Handle, double-ended

Kirsch Sinus Lift Instrument Set SMALL

for sinus floor elevation, modif. acc. to Dr. Kai Zwanzig

Responding to numerous customer requests, the proven Kirsch shapes were minimized during the further development of the new sinus lift instruments SMALL.

- The smaller instruments allow the preparation of a minimalized window.
- Regarding the angles of working tips, the previous Kirsch Instrument Set left nothing to be desired.

- The minimized analogous instruments now enable the practitioner to perform the preparation in a smaller area without injuring the Schneider's membrane.
- The **ZEPF**-Line-handle allows a safe, tactile guidance of the instruments.

Optionally, two instruments 41.868.05SL and 41.868.06SL with a longer 1st shaft are offered.



41.868.05S

Kirsch, Sinus 5, SMALL, Elevator 90°/90°/60°, Titanium **ZEPF**-Line Handle, double-ended



41.868.06S

Kirsch, Sinus 6, SMALL, Elevator 90°/90°, Titanium **ZEPF**-Line Handle, double-ended



41.868.07S

Kirsch, Sinus 7, SMALL, Plugger Ø 5.0 mm, Spoon 8 x 10 mm, Titanium **ZEPF**-Line Handle, double-ended



41.868.08S

Kirsch, Sinus 8, SMALL, convex Plugger Ø 4.5 mm, plane Plugger Ø 2.7 mm, Titanium **ZEPF**-Line Handle, double-ended



41.747.03S

Universal Raspatory, SMALL, exchangeable inserts



Kirsch Sinus Lift Instrument Set SMALL

for sinus floor elevation, modif. acc. to Dr. Kai Zwanzig

Optional Instrumentary



41.868.05SL

Kirsch, Sinus 5, SMALL, long shank,
Elevator 90°/90°/60°,
Titanium **ZEPF**-Line Handle, double-ended



41.868.06SL

Kirsch, Sinus 6, SMALL, long shank,
Elevator 90°/90°,
Titanium **ZEPF**-Line Handle, double-ended



AESTHETIC IS THE RESULT